

# HONEYWELL DUAL LAYER FACE COVER

## REUSABLE KNITTED FACE COVER WITH INSERT

This face mask product is designed to withstand multiple washes, and when used with a disposable insert provides the wearer with an extra layer of mouth and nose coverage beyond a common cloth face covering.

### PRODUCT FEATURES

- 4-way stretch minimizes gaps and enhances overall comfort
- Adjustable ear loops and nose clip
- Contoured design stands off the face to enhance breathability
- Hidden pocket designed for easy insert removal and replacement.
- Single-use inserts are disposable.
- Replace daily and dispose of inserts properly after use.
- Face mask maintains function and shape after multiple washings
- Remove insert prior to washing
- Stretchable 3D knitted fabric to fit multiple face profiles
- Ventilation through knitting pattern for improved breathability
- Colors: Dark Gray, Light Gray

### TECHNICAL PARAMETERS

- Face Cover with Protective Insert\*
  - ✓ Washable face cover for ease of reusable
  - ✓ Disposable inserts for single use



### SKU INFORMATION

PART NUMBER	DESCRIPTION
LM-500-WF	1 Dual-Layer Face Cover Dark Grey with Grey Noseclip + 12 Inserts
LM-500LG-WF	1 Dual-Layer Face Cover Light Grey with Grey Noseclip + 12 Inserts
LMF-500-12	12 Pack of Inserts

\*The insert has been shown by a third-party testing facility to block up to 97% of 3-micron-sized non-viable particles and 0.1-micron-sized aerosolized particles

#### For More Information

[www.safety.honeywell.com](http://www.safety.honeywell.com)

#### Honeywell PPE

9680 Old Bailes Rd.  
Fort Mill, SC 29707

WARNING: The product has not been FDA cleared or approved. The product has been authorized by FDA under an EUA for use as source control by the general public as well as by HCP in healthcare settings as to help prevent the spread of infection or illness during the COVID-19 pandemic. Source control refers to the use of a face mask over the mouth and nose to contain that individual's respiratory secretions to help prevent transmission from infected individuals who may or may not have symptoms of COVID-19. This product is authorized only for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of medical devices, including alternative products used as medical devices, during the COVID-19 outbreak, under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1) unless the authorization is terminated or revoked sooner. The face mask is intended to be used with single-use inserts. There has not been testing of the effectiveness of the disposable insert when placed in the mask. Accordingly, this product is not a respiratory protective device, nor a substitute for filtering face piece respirators or for surgical face masks. This product is not for medical use or for use in a high infection risk clinical setting. It is not intended to be used to provide liquid barrier protection or for high-risk aerosol generating procedures.

Update 2/2021

THE  
FUTURE  
IS  
WHAT  
WE  
MAKE IT

Honeywell