

This document is intended to help customers to select medical SPO2 sensors and adapters out of the Honeywell Healthcare Solutions GmbH portfolio that are compatible to the original sensors.

While Honeywell may provide information or engineering support for its products through Honeywell personnel, literature and website, it is the customer's sole responsibility to determine the suitability of the Honeywell product(s) for the customer's requirements. Specifications may change without notice. The information we supply is believed to be accurate and reliable as of this writing.

Check the Honeywell website for the latest version of this document and for technical details of the sensor specifications!

[Link to Cross Reference documents on website](https://hwll.co/envitec-cross-reference-list)

<https://hwll.co/envitec-cross-reference-list>

[Information about SPO2 sensors](https://automation.honeywell.com/de/de/products/sensing-solutions/sensors/pulse-oximetry-sensors)

<https://automation.honeywell.com/de/de/products/sensing-solutions/sensors/pulse-oximetry-sensors>

Contact: info@envitec.com

Appendix:

- Connector Guide
- SPO2Sensor compatibility list
- SPO2 Adapter list

	Revision:	Valid date:	Created:	Page – Seite:
006-35-XRL-Cross_Reference_List_SPO2-0.docx	0	2025-07-17	HB	1 / 1

CONNECTOR GUIDE

ENVITEC
by Honeywell



Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Armstrong Medical	AD-1000	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Armstrong Medical	AD-1000	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Armstrong Medical	AD-1000	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Armstrong Medical	AD-1000	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-1000	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Armstrong Medical	AD-2000	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Armstrong Medical	AD-2000	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Armstrong Medical	AD-2000	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Armstrong Medical	AD-2000	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-2000	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Armstrong Medical	AD-3000	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Armstrong Medical	AD-3000	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Armstrong Medical	AD-3000	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Armstrong Medical	AD-3000	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-3000	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Armstrong Medical	AD-4000	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Armstrong Medical	AD-4000	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Armstrong Medical	AD-4000	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Armstrong Medical	AD-4000	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-4000	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Armstrong Medical	AD-900	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Armstrong Medical	AD-900	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Armstrong Medical	AD-900	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Armstrong Medical	AD-900	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Armstrong Medical	AD-900	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Artema	Diascope	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	59831-0
Artema	Diascope	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	59831-0

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Artema	Diascope	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	59831-0
Artema	Diascope	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	59831-0
Baxter	Explorer/ Monarch	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Baxter	Explorer/ Monarch	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
BCI (Biochem Int'l)	3101	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3402	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3402	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3100	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3100	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3101	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3101	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3101	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	3180	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3180	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3180	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3300	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3300	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3301	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3301	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3302	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3302	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3303	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3303	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3304	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3304	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3304	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	3401	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	3401	not available	FingerClip	E06-00-1123	available	F-3222-30 BCI	3	C03	direct	
BCI (Biochem Int'l)	6100	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	6100	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	6100	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	8400 Capnocheck	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	9100	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	9100	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	9100	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Advisor	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Advisor	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Advisor	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Autocorr	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Autocorr	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Autocorr	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Mini Torr	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Mini Torr	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Mini Torr Plus	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	BCI® 3109
BCI (Biochem Int'l)	Spectro 10	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Spectro 10	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 10	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Spectro 20	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 20	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 1300	Disposable sensor	E08-00-0123	phased out	DA-2221-1 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
BCI (Biochem Int'l)	Spectro 30	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
BCI (Biochem Int'l)	Spectro 30	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Biosys Co., LTD	BPM-200	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Biosys Co., LTD	BPM-200	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Biosys Co., LTD	BPM-200	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Biosys Co., LTD	Vital Sign Monitor SENTRY	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
CAS Medical Systems Inc.	OscilloMate 9002		Disposable sensor	E08-00-0127	available	DA-2231-1 Nonin	0,45	C04	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
CAS Medical Systems Inc.	OscilloMate 9002	7000I	Disposable sensor	E08-00-0129	available	DI-2231-5 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	7000N	Disposable sensor	E08-00-0130	available	DN-2231-6 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	7000P	Disposable sensor	E08-00-0128	available	DP-2231-2 Nonin	0,45	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000 J	Wrap-Sensor	E09-00-0270	phased out	W-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000AA-1	FingerClip	E06-00-0216	available	F-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SL-1	SoftTip® large	E06-00-0765	available	R-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SL-1	SoftTip® medium	E06-00-1368	phased out	RM-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SL-1	SoftTip® plus large	E1001383	phased out	R-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SL-1	SoftTip® plus medium	E1001384	phased out	RM-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SS-1	SoftTip® small	E06-00-0814	available	RS-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9002	Nonin® 8000SS-1	SoftTip® plus small	E1001385	phased out	RS-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S		Disposable sensor	E08-00-0127	available	DA-2231-1 Nonin	0,45	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	7000I	Disposable sensor	E08-00-0129	available	DI-2231-5 Nonin	0,9	C04	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
CAS Medical Systems Inc.	OscilloMate 9302S	7000N	Disposable sensor	E08-00-0130	available	DN-2231-6 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	7000P	Disposable sensor	E08-00-0128	available	DP-2231-2 Nonin	0,45	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000 J	Wrap-Sensor	E09-00-0270	phased out	W-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000AA-1	FingerClip	E06-00-0216	available	F-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SL-1	SoftTip® large	E06-00-0765	available	R-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SL-1	SoftTip® medium	E06-00-1368	phased out	RM-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SL-1	SoftTip® plus large	E1001383	phased out	R-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SL-1	SoftTip® plus medium	E1001384	phased out	RM-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SS-1	SoftTip® small	E06-00-0814	available	RS-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9302S	Nonin® 8000SS-1	SoftTip® plus small	E1001385	phased out	RS-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303		Disposable sensor	E08-00-0127	available	DA-2231-1 Nonin	0,45	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	7000I	Disposable sensor	E08-00-0129	available	DI-2231-5 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	7000N	Disposable sensor	E08-00-0130	available	DN-2231-6 Nonin	0,9	C04	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
CAS Medical Systems Inc.	OscilloMate 9303	7000P	Disposable sensor	E08-00-0128	available	DP-2231-2 Nonin	0,45	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000 J	Wrap-Sensor	E09-00-0270	phased out	W-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000AA-1	FingerClip	E06-00-0216	available	F-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SL-1	SoftTip® large	E06-00-0765	available	R-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SL-1	SoftTip® medium	E06-00-1368	phased out	RM-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SL-1	SoftTip® plus large	E1001383	phased out	R-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SL-1	SoftTip® plus medium	E1001384	phased out	RM-3012plus Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SS-1	SoftTip® small	E06-00-0814	available	RS-3012 Nonin	0,9	C04	direct	
CAS Medical Systems Inc.	OscilloMate 9303	Nonin® 8000SS-1	SoftTip® plus small	E1001385	phased out	RS-3012plus Nonin	0,9	C04	direct	
Choice Electronic Technologie Co	300M	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	300M	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	300M	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	300M	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MD300I	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300I	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MD300K1	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MD300K1	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MD300K1	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M3	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-M7	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-S6	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Choice Electronic Technologie Co	MMED6000 DP-SF	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Colin	BP 88	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	use original cable
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP 88	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	use original cable

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-306	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-306	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	use original cable
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-508	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-508	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-808	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-88S	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	BP-88S	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	BP-88S	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Colin	Press-Mate SAT	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N.3000	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N.3000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.3000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N.390	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N.390	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-10	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	N-10	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-10	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	N-10	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	N-10	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	N-10	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-10	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-100 (w / High Sat)	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	E610-01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	E610-01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-100 (w / High Sat)	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-100 (w / High Sat)	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-180	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	N-180	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	N-180	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-180	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	N-180	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-180	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	N-190	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	N-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	N-190	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	N-190	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-190	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-20	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	N-20	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-20	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	N-20	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	N-20	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	N-20	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-20	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-20	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-200	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	E610-01
Covidien	N-200	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	E610-01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	E610-01
Covidien	N-200	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	E610-01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-200	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	E610-01
Covidien	N-200	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-200	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	E610-01
Covidien	N-30	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	N-30	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-30	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	N-30	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	N-30	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	N-30	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	N-30	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-6000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	N-6000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	N-6000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-190	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-190	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-190	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-190	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-190	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-190	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-195	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-195	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-195	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-195	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-195	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-290/5	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-290/5	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-290/5	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Covidien	NPB-3910	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-3910	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3910	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3910	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3910	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3910	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-3920	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3920	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3920	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3920	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3920	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3930	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-3930	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3930	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3930	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3930	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3930	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3930	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-3940	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-3940	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-3940	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-3940	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-3940	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-4000	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-4000	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-4000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-4000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-4000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-75	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Covidien	NPB-75	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® large	E06-00-0709	phased out	R-3212-31 Nellcor	3	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology) + EC8	FingerClip	E06-00-1075	phased out	F-3212-31 Nellcor	3,1	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus large	E1001374	available	R-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-75	Nellcor® DS100A (non OxiMax technology) + EC8	SoftTip® plus medium	E1001375	available	RM-3212-31plus Nellcor	3	C05	direct	
Covidien	NPB-75	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® D-YS (non OxiMax technology) + EC8	Y-Sensor	E09-00-0119	phased out	Y-2217-31 Nellcor	3,1	C05	direct	
Covidien	NPB-75	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Covidien	NPB-75	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Covidien	NPB-75	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Accutorr Plus (Nellcor® SpO2)	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
Datascope	Expert (Nellcor® SpO2)	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Expert (Nellcor® SpO2)	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Datascope	Expert (Nellcor® SpO2)	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Datascope	Passport (Nellcor® SpO2)	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Datascope	Passport (Nellcor® SpO2)	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Digicare	Digimax 5500	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	use original cable

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Digicare	Digimax 5500	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	indirect	use original cable
Digicare	Digimax 5500	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	indirect	use original cable
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	M35370

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Cicero	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	M35370
Dräger	Cicero	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	M35370
Dräger	Cicero	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Dräger Oxipac 2000	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Dräger Oxipac 2000	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2000	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 1301	Disposable sensor	E08-00-0124	phased out	DP-2221-2 BCI	0,45	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 1302	Disposable sensor	E08-00-0125	phased out	DI-2221-5 BCI	0,9	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 1303	Disposable sensor	E08-00-0126	phased out	DN-2221-6 BCI	0,9	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Dräger Oxipac 2500	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3078	Ear sensor	E06-00-0585	phased out	ES-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
Dräger	Dräger Oxipac 2500	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Dräger	Infinity® Delta	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® Delta	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Delta	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® DeltaXL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® Gamma	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® Gamma X XL	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Gamma X XL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® GammaXL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® GammaXL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor« D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® Kappa	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Kappa	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® KappaXL	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® KappaXL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® M300	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® M300	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433 (1m)
Dräger	Infinity® Vista	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® Vista	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® Vista	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	3368433

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Infinity® VistaXL	Nellcor® D-YS (non OxiMax technology)	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Infinity® VistaXL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	3368433
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	pre-amp

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Narkomed	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Narkomed	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	pre-amp
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	M35370

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Dräger Fingerclip Sensor MS13235	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	M35370
Dräger	Vitalert 1000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	M35370
EnviteC	ChipOx Modul		SoftTip® large	E06-00-0472	available	R-3227 EnviteC	1,2	C45	direct	
EnviteC	ChipOx Modul		FingerClip	E06-00-1268	available	F-3227 EnviteC	1,2	C45	direct	
EnviteC	ChipOx Modul		SoftTip® small	E06-00-1355	available	RS-3227 EnviteC	1,2	C45	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
EnviteC	ChipOx Modul		SoftTip® medium	E06-00-1371	available	RM-3227 EnviteC	1,2	C45	direct	
EnviteC	ChipOx Modul		Ear sensor	E06-00-1432	available	ES-3227 EnviteC	1,2	C45	direct	
EnviteC	ChipOx Modul		Y-Sensor	E08-00-0464	phased out	Y-3227 EnviteC	1,2	C45	direct	
EnviteC	ChipOx Modul		Wrap-Sensor	E09-00-0261	phased out	W-3227 EnviteC	1,2	C45	direct	
EnviteC	MySign® S		SoftTip® large	E1002015	available	R-3211-31 MySign®	2,95	C74	direct	
EnviteC	MySign® S		SoftTip® large	E1002016	available	R-3211-12 MySign®	1,2	C74	direct	
EnviteC	MySign® S		SoftTip® medium	E1002018	available	RM-3211-12 MySign®	1,2	C74	direct	
EnviteC	MySign® S		SoftTip® small	E1002020	phased out	RS-3211-12 MySign®	1,2	C74	direct	
EnviteC	MySign® S		SoftTip® plus large	E1002021	available	R-3211-31plus MySign®	2,95	C74	direct	
EnviteC	MySign® S		SoftTip® plus large	E1002022	available	R-3211-12plus MySign®	1,2	C74	direct	
EnviteC	MySign® S		FingerClip	E1002027	available	F-3211-31 MySign®	2,95	C74	direct	
EnviteC	MySign® S		FingerClip	E1002028	phased out	F-3211-12 MySign®	1,2	C74	direct	
EnviteC	MySign® S		Ear sensor	E1002031	phased out	ES-3211-31 MySign®	2,95	C74	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
EnviteC	MySign® S		Ear sensor	E1002032	phased out	ES-3211-12 MySign®	1,2	C74	direct	
EnviteC	MySign® S		SoftTip® large	E1002038	phased out	R-3212-9 MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		SoftTip® medium	E1002039	available	RM-3212-9 MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		SoftTip® small	E1002045	phased out	RS-3212-9 MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		FingerClip	E1002047	phased out	F-3212-9 MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		Disposable sensor	E1002049	phased out	DA-2211-1 MySign®	0,45	C05	indirect	X-4211-1
EnviteC	MySign® S		Ear sensor	E1002050	phased out	ES-3212-9 MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		Disposable sensor	E1002053	phased out	DP-2211-2 MySign®	0,45	C05	indirect	X-4211-1
EnviteC	MySign® S		SoftTip® plus large	E1002054	available	R-3212-9plus MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		SoftTip® plus medium	E1002055	available	RM-3212-9plus MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S		SoftTip® plus small	E1002056	available	RS-3212-9plus MySign®	0,9	C05	indirect	X-4211-1
EnviteC	MySign® S	BCI« 1303	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
EnviteC	OxiPen®		SoftTip® large	E06-00-0472	available	R-3227 EnviteC	1,2	C45	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
EnviteC	OxiPen®		FingerClip	E06-00-1268	available	F-3227 EnviteC	1,2	C45	direct	
EnviteC	OxiPen®		SoftTip® small	E06-00-1355	available	RS-3227 EnviteC	1,2	C45	direct	
EnviteC	OxiPen®		SoftTip® medium	E06-00-1371	available	RM-3227 EnviteC	1,2	C45	direct	
EnviteC	OxiPen®		Ear sensor	E06-00-1432	available	ES-3227 EnviteC	1,2	C45	direct	
EnviteC	OxiPen®		Y-Sensor	E08-00-0464	phased out	Y-3227 EnviteC	1,2	C45	direct	
EnviteC	OxiPen®		Wrap-Sensor	E09-00-0261	phased out	W-3227 EnviteC	1,2	C45	direct	
EnviteC	PulsiQuant 8040	BCI® 3043	Y-Sensor	E09-00-0086	available	Y-2227 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3044	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3044	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3044	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3044	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3444N	FingerClip	E06-00-0232	phased out	F-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3444N	SoftTip® small	E06-00-0583	phased out	RS-3222-12 BCI	1,2	C03	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
EnviteC	PulsiQuant 8040	BCI® 3444N	SoftTip® large	E06-00-0584	phased out	R-3222-12 BCI	1,2	C03	direct	
EnviteC	PulsiQuant 8040	BCI® 3444N	SoftTip® medium	E06-00-1370	phased out	RM-3222-12 BCI	1,2	C03	direct	
Fukuda Denshi	DS-5100E	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Fukuda Denshi	DS-5100E	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5100E	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® large	E06-00-0711	available	R-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E06-00-1079	available	ES-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
Fukuda Denshi	DS-5300W	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
Fukuda Denshi	DS-5300W	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
Fukuda Denshi	LX-5230		FingerClip	E06-00-1066	phased out	F-2182-3 Minolta	0,3	C45	direct	
GE Medical Systems	7000 Series	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems	7000 Series	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems	7000 Series	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	7000 Series	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	7000 Series	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	7000 Series	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	7000 Series	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	7000 Series	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	7000 Series	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	7000 Series	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	7000 Series	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	7000 Series	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	Corometrics 118	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 118	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 118	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 120	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 120	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 128	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 128	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 128F	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 128F	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 129	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 129F	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 129F	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Corometrics 510	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 510	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 511	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 511	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 511	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 556	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 570	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Corometrics 570	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	4033CAX
GE Medical Systems	Corometrics 570	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9700	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9700	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9710	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9710	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap 9720 Plus	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Compact	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap Compact	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap MPS	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap Select	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Critikon Dinamap Select	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	direct	
GE Medical Systems	Critikon Dinamap Select	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	direct	
GE Medical Systems	Dash	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dash	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dash	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Dinamap Pro 100	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 100	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 1000 V2	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 200	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 200	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 300	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 300	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	SCP-10 / MC-10

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Dinamap Pro 400	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Dinamap Pro 400	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	Eagle	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Eagle	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Eagle	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	Hellige Eagle 1000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle 1000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 1000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle 3000	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle SMK	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle SMK	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMK	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle SMU	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMU	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	303 443 58

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	Hellige Eagle SMV	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® D-Ys (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	Hellige Eagle SMV	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	303 443 58
GE Medical Systems	SOLAR	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	SOLAR	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	SOLAR	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 200SL	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 200SL	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 200SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 200SL	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 200SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 200SL	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 200SL	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 200SL	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 200SL	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM 200SL	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 200SL	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 250SL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 250SL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 250SL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 250SL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 250SL	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 400SL	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 400SL	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 400SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 400SL	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 400SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 400SL	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 400SL	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 400SL	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 400SL	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM 400SL	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 400SL	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 450SL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 450SL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 450SL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 450SL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 450SL	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 600SL	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 600SL	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 600SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 600SL	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 600SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 600SL	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 600SL	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 600SL	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 600SL	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM 600SL	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 600SL	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 650SL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 650SL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 650SL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 650SL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 650SL	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 800A	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 800A	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 800A	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800A	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800A	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800A	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800A	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800A	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 800A	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM 800A	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800A	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 800SL	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM 800SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800SL	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800SL	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800SL	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 800SL	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM 800SL	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM 800SL	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 800SL	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM 850A	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 850A	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 850A	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 850A	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 850A	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850A	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 850SL	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	SCP-10 / MC-10
GE Medical Systems	TRAM 850SL	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM 850SL	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM A.R.	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM A.R.	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	223 394 01
GE Medical Systems	TRAM A.R.	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM A.R.	Ohmeda® 6051-0000-035	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM A.R.	Ohmeda® 6051-0000-035	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM A.R.	Ohmeda® 6051-0000-035	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems	TRAM A.R.	Ohmeda® 6051-0000-035	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems	TRAM A.R.	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM A.R.	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems	TRAM A.R.	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM A.R.	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems	TRAM Modular	Nellcor® D-20	Disposable sensor	E08-00-0469	available	DP-2211-2 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® D-25	Disposable sensor	E08-00-0468	available	DA-2211-1 Nellcor	0,45	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	SoftTip® small	E06-00-0712	available	RS-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	FingerClip	E06-00-1074	available	F-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	SoftTip® medium	E06-00-1369	available	RM-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus large	E1001363	phased out	R-3212-9plus Nellcor	0,9	C05	indirect	223 394 01

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus medium	E1001366	phased out	RM-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS100A (non OxiMax technology)	SoftTip® plus small	E1001367	phased out	RS-3212-9plus Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® DS-100A (non OxiMax technology)	SoftTip® large	E1001975	available	R-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® D-YS (non OxiMax technology)	Y-Sensor	E09-00-0087	phased out	Y-2217-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® D-YSE (non OxiMax technology)	Ear sensor	E1001974	available	ES-3217-9 GE Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® I-20	Disposable sensor	E08-00-0470	available	DI-2211-5 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® N-25	Disposable sensor	E08-00-0471	available	DN-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® N-25	Disposable sensor	E1000947	available	DW-2211-6 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems	TRAM Modular	Nellcor® OXI-A/N Oxiband	Wrap-Sensor	E09-00-0266	phased out	W-3212-9 Nellcor	0,9	C05	indirect	223 394 01
GE Medical Systems (Datex-Ohmeda)	AS/3	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	AS/3	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	AS/3	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Capnomac	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Capnomac	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap 5	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap I	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap II	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Cardiicap/5 with N-XOSAT option S/5	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CG/2S	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CG/2S	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CG/2S	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CG/2S	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CG/2S	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CG/2S	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH/S	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH/S	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CH/S	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CH/S	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH/S	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH-RS	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH-RS	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH-RS	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH-RS	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CH-RS	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CS/3	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CS/3	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	SoftTip® medium	E1001246	phased out	RM-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	SoftTip® large	E1001247	phased out	R-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	FingerClip	E1001248	phased out	F-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	SoftTip® small	E1001545	phased out	RS-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-F-UN + OXY-ES3	SoftTip® large	E1001819	available	R-3513-40 GE Datex-Ohmeda	4	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Light Monitor	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Light Monitor	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Light Monitor	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Light Monitor	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Light Monitor	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Light Monitor	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	M-Estpr	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	M-Estpr	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	M-Estpr	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	M-Estpr	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	M-Estpr	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Module E-PREST	OXY-F-UN + OXY-ES3	SoftTip® medium	E1001246	phased out	RM-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PREST	OXY-F-UN + OXY-ES3	SoftTip® large	E1001247	phased out	R-3513-30 GE Datex-Ohmeda	3	C64	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Module E-PREST	OXY-F-UN + OXY-ES3	FingerClip	E1001248	phased out	F-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PREST	OXY-F-UN + OXY-ES3	SoftTip® small	E1001545	phased out	RS-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PREST	OXY-F-UN + OXY-ES3	SoftTip® large	E1001819	available	R-3513-40 GE Datex-Ohmeda	4	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSM	OXY-F-UN + OXY-ES3	SoftTip® medium	E1001246	phased out	RM-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSM	OXY-F-UN + OXY-ES3	SoftTip® large	E1001247	phased out	R-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSM	OXY-F-UN + OXY-ES3	FingerClip	E1001248	phased out	F-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSM	OXY-F-UN + OXY-ES3	SoftTip® small	E1001545	phased out	RS-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSM	OXY-F-UN + OXY-ES3	SoftTip® large	E1001819	available	R-3513-40 GE Datex-Ohmeda	4	C64	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Module E-PSMP	OXY-F-UN + OXY-ES3	SoftTip® medium	E1001246	phased out	RM-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSMP	OXY-F-UN + OXY-ES3	SoftTip® large	E1001247	phased out	R-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSMP	OXY-F-UN + OXY-ES3	FingerClip	E1001248	phased out	F-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSMP	OXY-F-UN + OXY-ES3	SoftTip® small	E1001545	phased out	RS-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-PSMP	OXY-F-UN + OXY-ES3	SoftTip® large	E1001819	available	R-3513-40 GE Datex-Ohmeda	4	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-REST	OXY-F-UN + OXY-ES3	SoftTip® medium	E1001246	phased out	RM-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-REST	OXY-F-UN + OXY-ES3	SoftTip® large	E1001247	phased out	R-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-REST	OXY-F-UN + OXY-ES3	FingerClip	E1001248	phased out	F-3513-30 GE Datex-Ohmeda	3	C64	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Module E-REST	OXY-F-UN + OXY-ES3	SoftTip® small	E1001545	phased out	RS-3513-30 GE Datex-Ohmeda	3	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Module E-REST	OXY-F-UN + OXY-ES3	SoftTip® large	E1001819	available	R-3513-40 GE Datex-Ohmeda	4	C64	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	not available	SoftTip® small	E06-00-0817	available	RS-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	not available	SoftTip® medium	E06-00-0824	available	RM-3412-9 Ohmeda	0,9	C01	indirect	X-2241-8
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F1-H	SoftTip® small	E06-00-0822	available	RS-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F1-H	FingerClip	E06-00-1080	phased out	F-3412-10 Ohmeda	1	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	FingerClip	E06-00-0264	phased out	F-3412 Ohmeda	3,6	C22	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	SoftTip® large	E06-00-0715	available	R-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	SoftTip® small	E06-00-0716	phased out	RS-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	SoftTip® large	E06-00-0771	phased out	R-3412-10 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	SoftTip® medium	E06-00-1377	phased out	RM-3412 Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F4-H	SoftTip® plus small	E1001394	available	RS-3412plus Ohmeda	3	C22	direct	
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4005

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4005
GE Medical Systems (Datex-Ohmeda)	Oscar II	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Oscar II	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Oscar II	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Oscar II	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Oscar II	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Oscar II	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	OSP-200	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	OSP-200	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	OSP-200	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	OSP-200	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	OSP-200	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	S/5 Light	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	S/5 Light	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	S/5 Light	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	S/5 Light	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	S/5 Light	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	Ultima	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Ultima	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Ultima	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Ultima	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Ultima	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTS00	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTS00	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTS00	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSIO	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSIO	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSIO	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSIO	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSIO	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSV I	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSV I	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSV I	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSV I	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSV I	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	Datex® SAS-F	FingerClip	E06-00-0268	phased out	F-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSVO	Datex® SAS-F	SoftTip® large	E06-00-0766	available	R-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSVO	Datex® SAS-F	SoftTip® medium	E06-00-1378	phased out	RM-3512-9 Datex	0,9	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSVO	not available	SoftTip® small	E06-00-0729	phased out	RS-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSVO	not available	SoftTip® small	E06-00-0815	phased out	RS-3512-9 Datex	0,9	C01	indirect	OXY-C3

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-AF-10	Disposable sensor	E1001712	available	DW-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-AF-10	Disposable sensor	E1001802	phased out	DN-3003-6 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-AF-10	Disposable sensor	E1001803	available	DP-3003-2 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-AF-10	Disposable sensor	E1001804	phased out	DI-3003-5 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-AP-10	Disposable sensor	E1001679	phased out	DA-3003-1 FlexMon	0,45	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-E-UN	Ear sensor	E1001674	phased out	ES-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F4-N	FingerClip	E06-00-0270	phased out	F-3512-40 Datex	3,6	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F4-N	SoftTip® large	E06-00-0727	phased out	R-3512-40 Datex	3	C14	direct	

Original manufacturer	Original model	Original article reference	Sensor type	Honeywell part no.	Product status	Honeywell sensor	Cable length [m]	plug	Connection option	Required original accessories
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F4-N	SoftTip® medium	E06-00-1380	phased out	RM-3512-40 Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F4-N	SoftTip® plus large	E1001377	available	R-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F4-N	SoftTip® plus medium	E1001378	phased out	RM-3512-40plus Datex	3	C14	direct	
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F-DB	FingerClip	E06-00-1119	phased out	F-3512-20 Datex	2	C01	indirect	OXY-C3
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-F-UN	FingerClip	E1001673	phased out	F-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	TS-SA-D	SoftTip® large	E1001709	phased out	R-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	TS-SA-D	SoftTip® medium	E1001710	phased out	RM-3003-9 FlexMon	0,9	C69	indirect	X-4004
GE Medical Systems (Datex-Ohmeda)	ULTSVO	TS-SA-D	SoftTip® small	E1001711	phased out	RS-3003-9 FlexMon	0,9	C69	indirect	X-4004

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
GE Medical Systems (Datex-Ohmeda)	Cardicap/5 with N-XOSAT option S/5	OXY-ES-3	FlexMon	E1001676	phased out	X-4003 GE Datex Ohmeda	3	C64	C70
GE Medical Systems (Datex-Ohmeda)	Datex-Ohmeda S/5 FM port	OXY-ES-3	FlexMon	E1001676	phased out	X-4003 GE Datex Ohmeda	3	C64	C70
GE Medical Systems (Datex-Ohmeda)	AS/3	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	CG/2S	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	CH-RS	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	CH/S	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	CS/3	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Capnomac Ultima	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Capnomac	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Cardicap 5	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Cardicap I	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Cardicap II	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Light Monitor	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
GE Medical Systems (Datex-Ohmeda)	M-Estpr	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	OSP-200	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Oscar II	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	S/5 Light	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus II	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Satelite Trans Plus	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	ULTS00	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	ULTSIO	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	ULTSV I	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	ULTSVO	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	Ultima	OXY-SL-3	FlexMon	E1001729	phased out	X-4004 Datex	3	C14	C70
GE Medical Systems (Datex-Ohmeda)	AS/3 modular with M-OSAT module	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
GE Medical Systems (Datex-Ohmeda)	CS/3 modular with M-OSAT module	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3700	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3710	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3760	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3770	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3775	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3800	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 3900P	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 4700	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda 5250RGM	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda Biox 3740	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda OxiCap	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
GE Medical Systems (Datex-Ohmeda)	Ohmeda Rascal II	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems (Datex-Ohmeda)	Ohmeda TuffSat	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	7000 Series	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 200SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 250SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 400SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 450SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 600SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 650SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 800A	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 800SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 850A	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM 850SL	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM A.R.	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
GE Medical Systems	TRAM Modular	OXY-OL3	FlexMon	E1001730	available	X-4005 Ohmeda	3	C22	C70
CAS Medical Systems Inc.	OscilloMate 9002	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
CAS Medical Systems Inc.	OscilloMate 9002	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
CAS Medical Systems Inc.	OscilloMate 9302S	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
CAS Medical Systems Inc.	OscilloMate 9302S	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
CAS Medical Systems Inc.	OscilloMate 9303	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
CAS Medical Systems Inc.	OscilloMate 9303	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Healthdyne	930	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Healthdyne	930	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Healthdyne	9705	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Healthdyne	9705	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Healthdyne	970S	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Healthdyne	970S	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Healthdyne	970SE	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Healthdyne	970SE	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Innomed	Innocare-T0	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Innomed	Innocare-T0	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Innomed	Oxycard	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Innomed	Oxycard	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Nonin	7500	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	7500	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	8500	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	8500	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	9840	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	9840	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	Avant® 2120	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	Avant® 2120	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	Avant® 9600	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	Avant® 9600	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	Avant® 9700	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60
Nonin	Avant® 9700	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Nonin	PamSAT® 2500	8500I	Extension and adapter cable	E06-00-0350	phased out	X-4012 Nonin	1	C04	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Nonin	PamSAT® 2500	8500VI	Extension and adapter cable	E06-00-0351	available	X-4012-1 Nonin	3	C04	C60
Baxter	Explorer/Monarch	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Biosys Co., LTD	BPM-200	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Biosys Co., LTD	Vital Sign Monitor SENTRY	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	N-10	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	N-180	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	N-190	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	N-20	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	N-30	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-190	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-195	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-3910	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-3920	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Covidien	NPB-3930	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-3940	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-4000	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-75	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Datascope	Expert (Nellcor® SpO2)	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Datascope	Passport (Nellcor® SpO2)	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Fukuda Denshi	DS-5100E	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Fukuda Denshi	DS-5300W	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
GE Medical Systems	Critikon Dinamap 9700	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
GE Medical Systems	Critikon Dinamap 9710	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
GE Medical Systems	Critikon Dinamap 9720 Plus	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
GE Medical Systems	Critikon Dinamap Compact	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
GE Medical Systems	Critikon Dinamap MPS	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
GE Medical Systems	Critikon Dinamap Select	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Ivy Biomedical	405P	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Ivy Biomedical	405T	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
M.D.E.	Escort II OPT30	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
M.D.E.	Prism	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Medsource	BPM-200	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Mennen	Horizon XL (Nellcor® SpO2)	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Mennen	Mercury	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Nihon Kohden	AL-900P Module	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Nihon Kohden	BSM-1102	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Nihon Kohden	BSM-9510	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Nihon Kohden	OGS-2000	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Welch Allyn Protocol	Propaq 106-EL	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Welch Allyn Protocol	Propaq 202-EL	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Welch Allyn Protocol	Propaq 204-EL	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Welch Allyn Protocol	Propaq 206-EL	EC-8	Extension and adapter cable	E06-00-0353	phased out	X-4212 Nellcor	2,4	C05	C60
Covidien	NPB-40 (non-Oximax)	EC-4	Extension and adapter cable	E06-00-0354	phased out	X-4212-1 Nellcor	1,2	C05	C60
Colin	BP 88	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Colin	BP-306	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Colin	BP-508	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Colin	BP-808	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Colin	BP-88S	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Colin	Press-Mate SAT	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Covidien	N-6000	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Covidien	N.3000	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Covidien	N.390	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Covidien	NPB-290/5	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Datascope	Accutorr Plus (Nellcor® SpO2)	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
GE Medical Systems	Dinamap Pro 100	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
GE Medical Systems	Dinamap Pro 1000 V2	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
GE Medical Systems	Dinamap Pro 200	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
GE Medical Systems	Dinamap Pro 300	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
GE Medical Systems	Dinamap Pro 400	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Nihon Kohden	BSM-4104A	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Nihon Kohden	BSM-4114A	SCP-10 / MC-10	Extension and adapter cable	E06-00-0355	phased out	X-4213 Nellcor	3	C41	C60
Philips Medical Systems	78352C	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	78354C	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	78834C	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Philips Medical Systems	M1020A	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M1025A	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M1025B	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M1350B	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M1722A/B	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M1732A/B	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M2475B	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	ViridaA1	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	ViridaA3	M1940A	Extension and adapter cable	E06-00-1147	phased out	X-5414 HP/Philips	2	C17	C62
Philips Medical Systems	M4529A	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62
Philips Medical Systems	MP20	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62
Philips Medical Systems	MP30	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62
Philips Medical Systems	MP60	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62

Original manufacturer	Original model	Original article reference	Extension / adapter cable	Honeywell part number	Product status	Honeywell adapter cable	Cable length [m]	Plug	Socket
Philips Medical Systems	MP70	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62
Philips Medical Systems	MP90	M1941A	Extension and adapter cable	E06-00-1194	available	X-5415 HP/Philips	2	C19	C62
Philips Medical Systems	M4529A	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
Philips Medical Systems	MP20	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
Philips Medical Systems	MP30	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
Philips Medical Systems	MP60	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
Philips Medical Systems	MP70	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
Philips Medical Systems	MP90	M1943A	Extension and adapter cable	E06-00-1196	phased out	X-5416 HP/Philips	1,1	C19	C60
EnviteC	ChipOx Modul		Extension and adapter cable	E06-00-1198	phased out	X-3227-12 EnviteC	1,2	C45	C63
EnviteC	OxiPen®		Extension and adapter cable	E06-00-1198	phased out	X-3227-12 EnviteC	1,2	C45	C63
EnviteC	MySign® S		Extension and adapter cable	E1002030	phased out	X-4211-1 MySign®	1	C74	C60